## English

GENERAL CLINICAL CHEMISTRY PROGRAMME: RQ9112/a, RQ9112/b, RQ9112/c

### CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to [www.riqas.net](http://www.riqas.net) to confirm the exact date on which you received this kit.

### CHARACTERISTICS

The pack contains lyophilised 5 ml samples. The vials are labelled with the sample number.

### PREPARATION/STORAGE/STABILITY

The samples are sealed under vacuum. Open the vial very carefully, avoiding any loss of material and using a **calibrated pipette**, reconstitute with an accurately measured **5 ml volume** of freshly double distilled water at **+20°C to +25°C**. Replace the rubber stopper, close vial and leave to stand for **60 minutes out of bright light** before use. Ensure that the lyophilised contents are completely dissolved by swirling gently. Do not shake the vial. Please analyse the sample within 2 days of reconstitution. The samples should be treated in the same way as patient samples. The samples should be stored at 2 - 8°C when not in use.

**Prostatic and Total Acid Phosphatase:** The samples should be stabilised by adding 1 drop of acetic acid solution to 1 ml of serum.

**Bilirubin** in the serum is light sensitive so the serum should be **stored in the dark**.

**NEFA** is stable for 24 hours at +2°C to +8°C following reconstitution.

**EGFR:** The scenarios required for calculation of this parameter are available on RIQASNet – RIQAS Kit Inserts.

Bacterial contamination of the reconstituted serum will cause reductions in the stability of many components.

### SAFETY

The serum is human based. It has been source tested and found to be negative for HBsAg and antibodies to HIV. For complete protection however, it is recommended that the serum be handled as carefully as patient specimens.

For **IN VITRO** use only.

**\* IMPORTANT NOTE:** Results must arrive at RIQAS by **17:00 HRS GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. Late results will not be accepted after the final date for the next sample.

GENERAL CLINICAL CHEMISTRY PROGRAMME / PROGRAMME CHIMIE CLINIQUE / PROGRAMMA PER CHIMICA CLINICA GENERALE /  
 PROGRAMA DE QUIMICA CLINICA GENERAL / PROGRAMA DE QUÍMICA CLÍNICA / ALMENT KLINISK KEMIPROGRAM /  
 ALLMÄNT KLINISKT KEMISKT PROGRAM / ΓΕΝΙΚΟ ΠΡΟΓΡΑΜΜΑ ΚΛΙΝΙΚΗΣ ΧΗΜΕΙΑΣ /  
 GENERALNY PROGRAM CHEMII KLINICZNEJ / 通用临床化学分析程序 / GENEL BİYOKİMYA PROGRAMI

RETURN OF RESULTS / DATES de RETOUR / INVIO DEI RISULTATI / ENVIO DE RESULTADOS / PLAN FOR INSENDELSE AF / CYKEL FÖR /  
 ΚΥΚΛΟΣ ΕΠΙΣΤΡΟΦΗΣ / PRZESYŁANIE WYNIKÓW / 回復結果 / SONUÇ GÖNDERİM TARİHLERİ / ODOSLANIE VÝSLEDKOV

CYCLE 68 / CICLO 68 / RESULTATER 68 / RESULTATRETUR 68 /  
 ΑΠΟΤΕΛΕΣΜΑΤΩΝ 68 / CYKL ZWROTU WYNIKÓW 68 / 返回結果周期 68 / DÖNEM 68 / CYKLUS 68

SAMPLE NO / ECHANTILLON NO. / N. CAMPIONE / N° DE MUESTRA / AMOSTRA Nº. / PRØVENR. / ΑΡ. ΔΕΙΓΜΑΤΟΣ / PROV NR. / ΗΜΕΡΟΜΗΝΙΑ/ NR PRÓBKÍ / 样品号 / ÖRNEK NO / Č. VZORKY	RECOMMENDED ANALYSIS DATE/ DATE RECOMMANDE POUR L' ANALYSE/ DATA DI ANALISI RACCOMANDATA/ FECHA DE ANALISIS RECOMENDADA/ DATA RECOMENDADA PARA ANÁLISE/ ANBEFALET ANALYSEDATO/ REKOMMENDERAT ANALYSDATUM/ ΣΥΝΙΣΤΩΜΕΝΗ ΗΜΕΡΟΜΗΝΙΑ ΑΝΑΛΥΣΗΣ ZALECANA DATA WYKONANIA OZNACZENIA/ 推荐分析日期/ ÖNERİLEN ANALİZ TARİHİ/ ODPORÚČANÝ DÁTUM ANALÝZY	* FINAL DATE/ *DATE LIMITE/ * DATA FINALE/ *FECHA FINAL/ * DATA FINAL/ * SLUTDATO/ * SLUTDATUM/ *TEΛΙΚΗ * DATA KOŃCOWA/ *最终日期/ *SON SONUÇ GÖNDERİM TARİHİ/ *FINÁLNY DÁTUM
---	---	---

1	13.03.23	20.03.23
2	27.03.23	03.04.23
3	10.04.23	17.04.23
4	24.04.23	01.05.23
5	08.05.23	15.05.23
6	22.05.23	29.05.23
7	05.06.23	12.06.23
8	19.06.23	26.06.23
9	03.07.23	10.07.23
10	17.07.23	24.07.23
11	31.07.23	07.08.23
12	14.08.23	21.08.23
13	28.08.23	04.09.23

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