

URINARY PROTEIN STANDARD (UP CAL)

LOT NUMBER: 291UP

SIZE: 1 x 5.5 ml

PREPARATION

Contents ready for use. Stable up to expiry quoted on the vial.

CONCENTRATION

0.98 g/l

SAFETY PRECAUTIONS

From Safety Data Sheet.

EC REP Randox Teoranta, Meenmore,
Dungloe, Donegal,
F94 TV06, Ireland

08 Sep '20 ne

THIS PAGE IS INTENTIONALLY BLANK