

## QP-H :Mean & expected ranges

**LOT** 22151



2024-10

Applicable to FDC analyzer using slides.  
 Applicable à l'analyseur FDC utilisant des lames.  
 Applicabile ad analizzatori FDC mediante slide.  
 Anwendbar auf den FDC-Analysator mit Objektträgern.  
 Aplicable al analizador FDC con slides.  
 Aplicável ao analisador FDC utilizando lâminas.  
 Применимо к анализатору FDC на основе слайдов.  
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.5	15.8 ~ 19.2	mmol/L
	315	285 ~ 345	mg/dL
BUN	18.31	16.53 ~ 20.10	mmol/L
	51.3	46.3 ~ 56.3	mg/dL
UA	625	571 ~ 672	µmol/L
	10.5	9.6 ~ 11.3	mg/dL
TCHO	4.27	3.80 ~ 4.73	mmol/L
	165	147 ~ 183	mg/dL
TG	2.52	2.25 ~ 2.80	mmol/L
	223	199 ~ 248	mg/dL
CRE	476	423 ~ 530	µmol/L
	5.39	4.78 ~ 6.00	mg/dL
TP	72	66 ~ 79	g/L
	7.2	6.6 ~ 7.9	g/dL
ALB	47	40 ~ 53	g/L
	4.7	4.0 ~ 5.3	g/dL

TBIL	77	62 ~ 91	$\mu\text{mol/L}$
	4.5	3.6 ~ 5.3	mg/dL
HDL-C	...	...	...
	...	...	...
Ca	3.10	2.70 ~ 3.48	mmol/L
	12.4	10.8 ~ 13.9	mg/dL
IP	2.00	1.78 ~ 2.20	mmol/L
	6.2	5.5 ~ 6.8	mg/dL
DBIL	...	...	...
	...	...	...
Mg	...	...	...
	...	...	...
GGT	709	605 ~ 813	U/L
GOT/AST	94	83 ~ 104	U/L
GPT/ALT	104	92 ~ 117	U/L
CPK	285	247 ~ 324	U/L
LDH	420	362 ~ 477	U/L
ALP	604	521 ~ 687	U/L
AMYL	202	174 ~ 229	U/L
LAP	76	61 ~ 91	U/L
CHE	346	300 ~ 391	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<http://www.fujifilm.com/products/medical/>



**FUJIFILM Corporation**

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN