

QP-H :Mean & expected ranges

LOT 22152



2024.12

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.5	15.8 ~ 19.2	mmol/L
	315	285 ~ 345	mg/dL
BUN	18.35	16.56 ~ 20.13	mmol/L
	51.4	46.4 ~ 56.4	mg/dL
UA	625	571 ~ 672	µmol/L
	10.5	9.6 ~ 11.3	mg/dL
TCHO	4.32	3.85 ~ 4.78	mmol/L
	167	149 ~ 185	mg/dL
TG	2.52	2.25 ~ 2.79	mmol/L
	223	199 ~ 247	mg/dL
CRE	479	424 ~ 534	µmol/L
	5.42	4.80 ~ 6.04	mg/dL
TP	72	66 ~ 79	g/L
	7.2	6.6 ~ 7.9	g/dL
ALB	47	41 ~ 54	g/L
	4.7	4.1 ~ 5.4	g/dL

TBIL	68	55 ~ 82	μmol/L
	4.0	3.2 ~ 4.8	mg/dL
HDL-C

Ca	3.10	2.70 ~ 3.48	mmol/L
	12.4	10.8 ~ 13.9	mg/dL
IP	2.00	1.81 ~ 2.23	mmol/L
	6.2	5.6 ~ 6.9	mg/dL
DBIL

Mg

GGT	700	597 ~ 803	U/L
GOT/AST	95	85 ~ 106	U/L
GPT/ALT	106	94 ~ 119	U/L
CPK	283	244 ~ 321	U/L
LDH	476	412 ~ 540	U/L
ALP	615	531 ~ 699	U/L
AMYL	205	177 ~ 232	U/L
LAP	76	60 ~ 91	U/L
CHE	346	301 ~ 392	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN