

QP-H :Mean & expected ranges

LOT 22254


2026.04

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.7	16.0 ~ 19.4	mmol/L
	319	288 ~ 349	mg/dL
BUN	18.31	16.53 ~ 20.10	mmol/L
	51.3	46.3 ~ 56.3	mg/dL
UA	619	571 ~ 672	µmol/L
	10.4	9.6 ~ 11.3	mg/dL
TCHO	4.24	3.78 ~ 4.71	mmol/L
	164	146 ~ 182	mg/dL
TG	2.50	2.22 ~ 2.78	mmol/L
	221	197 ~ 246	mg/dL
CRE	480	425 ~ 535	µmol/L
	5.43	4.81 ~ 6.05	mg/dL
TP	72	66 ~ 78	g/L
	7.2	6.6 ~ 7.8	g/dL
ALB	46	40 ~ 53	g/L
	4.6	4.0 ~ 5.3	g/dL

TBIL	70	56 ~ 86	μmol/L
	4.1	3.3 ~ 5.0	mg/dL
HDL-C

Ca	3.08	2.70 ~ 3.45	mmol/L
	12.3	10.8 ~ 13.8	mg/dL
IP	2.00	1.78 ~ 2.20	mmol/L
	6.2	5.5 ~ 6.8	mg/dL
DBIL

Mg

GGT	711	607 ~ 815	U/L
GOT/AST	90	80 ~ 100	U/L
GPT/ALT	101	89 ~ 112	U/L
CPK	268	232 ~ 304	U/L
LDH	424	366 ~ 481	U/L
ALP	626	541 ~ 712	U/L
AMYL	205	177 ~ 232	U/L
LAP	74	59 ~ 89	U/L
CHE	347	301 ~ 393	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN