

QP-H :Mean & expected ranges

LOT 22355

 2026.08

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicável ao analisador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.9	16.2 ~ 19.6	mmol/L
	322	291 ~ 353	mg/dL
BUN	17.74	15.99 ~ 19.49	mmol/L
	49.7	44.8 ~ 54.6	mg/dL
UA	648	595 ~ 702	µmol/L
	10.9	10.0 ~ 11.8	mg/dL
TCHO	4.16	3.70 ~ 4.63	mmol/L
	161	143 ~ 179	mg/dL
TG	2.64	2.36 ~ 2.92	mmol/L
	234	209 ~ 259	mg/dL
CRE	452	400 ~ 504	µmol/L
	5.11	4.52 ~ 5.70	mg/dL
TP	73	67 ~ 79	g/L
	7.3	6.7 ~ 7.9	g/dL
ALB	46	40 ~ 53	g/L
	4.6	4.0 ~ 5.3	g/dL

TBIL	72	58 ~ 86	$\mu\text{mol/L}$
	4.2	3.4 ~ 5.0	mg/dL
HDL-C

Ca	3.40	3.03 ~ 3.80	mmol/L
	13.6	12.1 ~ 15.2	mg/dL
IP	2.07	1.84 ~ 2.29	mmol/L
	6.4	5.7 ~ 7.1	mg/dL
DBIL

Mg

GGT	683	582 ~ 784	U/L
GOT/AST	106	94 ~ 118	U/L
GPT/ALT	115	102 ~ 129	U/L
CPK	272	235 ~ 308	U/L
LDH	387	333 ~ 440	U/L
ALP	635	548 ~ 721	U/L
AMYL	172	147 ~ 197	U/L
LAP	81	66 ~ 96	U/L
CHE	337	292 ~ 381	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları" nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN