

QP-L :Mean & expected ranges

LOT 21151

 2024-08

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicável ao analisador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item	Mean	Expected range	Unit
Élément de test	Moyenne	Plage attendue	Unité
Campione	Valore medio	Intervallo previsto	Unità
Testelement	Mittelwert	Erwarteter Bereich	Einheit
Elemento de prueba	Media	Rango esperado	Unidad
Item de teste	Média	Intervalo esperado	Unidade
Объект испытания	Среднее	Ожидаемые пределы	Единица
Test ögesi	Ortalama	Beklenen aralık	Birim
GLU	5.9	5.2 ~ 6.6	mmol/L
	106	94 ~ 118	mg/dL
BUN	5.11	4.43 ~ 5.78	mmol/L
	14.3	12.4 ~ 16.2	mg/dL
UA	202	172 ~ 232	µmol/L
	3.4	2.9 ~ 3.9	mg/dL
TCHO	3.52	3.13 ~ 3.90	mmol/L
	136	121 ~ 151	mg/dL
TG	1.42	1.21 ~ 1.64	mmol/L
	126	107 ~ 145	mg/dL
CRE	88	54 ~ 123	µmol/L
	1.00	0.61 ~ 1.39	mg/dL
TP	59	53 ~ 65	g/L
	5.9	5.3 ~ 6.5	g/dL
ALB	36	30 ~ 42	g/L
	3.6	3.0 ~ 4.2	g/dL

TBIL	12	7 ~ 17	$\mu\text{mol/L}$
	0.7	0.4 ~ 1.0	mg/dL
HDL-C	1.09	0.88 ~ 1.29	mmol/L
	42	34 ~ 50	mg/dL
Ca	2.18	1.85 ~ 2.50	mmol/L
	8.7	7.4 ~ 10.0	mg/dL
IP	1.13	0.97 ~ 1.26	mmol/L
	3.5	3.0 ~ 3.9	mg/dL
DBIL	6.8	1.7 ~ 10.3	$\mu\text{mol/L}$
	0.4	0.1 ~ 0.6	mg/dL
Mg	0.66	0.53 ~ 0.78	mmol/L
	1.6	1.3 ~ 1.9	mg/dL
GGT	63	52 ~ 74	U/L
GOT/AST	36	30 ~ 42	U/L
GPT/ALT	39	33 ~ 45	U/L
CPK	64	48 ~ 80	U/L
LDH	180	151 ~ 209	U/L
ALP	186	154 ~ 219	U/L
AMYL	62	50 ~ 74	U/L
LAP	48	33 ~ 63	U/L
CHE	283	246 ~ 321	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<http://www.fujifilm.com/products/medical/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN