

QP-L :Mean & expected ranges

LOT 21152

 2024-12

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	6.0	5.3 ~ 6.7	mmol/L
	108	96 ~ 120	mg/dL
BUN	5.14	4.46 ~ 5.82	mmol/L
	14.4	12.5 ~ 16.3	mg/dL
UA	202	172 ~ 232	µmol/L
	3.4	2.9 ~ 3.9	mg/dL
TCHO	3.54	3.13 ~ 3.93	mmol/L
	137	121 ~ 152	mg/dL
TG	1.46	1.23 ~ 1.68	mmol/L
	129	109 ~ 149	mg/dL
CRE	88	54 ~ 123	µmol/L
	1.00	0.61 ~ 1.39	mg/dL
TP	59	53 ~ 65	g/L
	5.9	5.3 ~ 6.5	g/dL
ALB	36	30 ~ 42	g/L
	3.6	3.0 ~ 4.2	g/dL

TBIL	14	9 ~ 19	μmol/L
	0.8	0.5 ~ 1.1	mg/dL
HDL-C	1.09	0.88 ~ 1.29	mmol/L
	42	34 ~ 50	mg/dL
Ca	2.18	1.83 ~ 2.50	mmol/L
	8.7	7.3 ~ 10.0	mg/dL
IP	1.16	1.00 ~ 1.29	mmol/L
	3.6	3.1 ~ 4.0	mg/dL
DBIL	8.6	3.4 ~ 12.0	μmol/L
	0.5	0.2 ~ 0.7	mg/dL
Mg	0.66	0.53 ~ 0.78	mmol/L
	1.6	1.3 ~ 1.9	mg/dL
GGT	64	53 ~ 75	U/L
GOT/AST	32	26 ~ 38	U/L
GPT/ALT	38	32 ~ 45	U/L
CPK	71	55 ~ 87	U/L
LDH	168	141 ~ 196	U/L
ALP	186	154 ~ 218	U/L
AMYL	61	49 ~ 73	U/L
LAP	49	34 ~ 64	U/L
CHE	282	244 ~ 319	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları" nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN