

## QP-L :Mean & expected ranges

**LOT** 21254


2026-01

Applicable to FDC analyzer using slides.  
 Applicable à l'analyseur FDC utilisant des lames.  
 Applicabile ad analizzatori FDC mediante slide.  
 Anwendbar auf den FDC-Analysator mit Objektträgern.  
 Aplicable al analizador FDC con slides.  
 Aplicável ao analisador FDC utilizando lâminas.  
 Применимо к анализатору FDC на основе слайдов.  
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	5.9	5.3 ~ 6.6	mmol/L
	107	95 ~ 119	mg/dL
BUN	5.14	4.46 ~ 5.82	mmol/L
	14.4	12.5 ~ 16.3	mg/dL
UA	208	172 ~ 238	µmol/L
	3.5	2.9 ~ 4.0	mg/dL
TCHO	3.47	3.08 ~ 3.85	mmol/L
	134	119 ~ 149	mg/dL
TG	1.42	1.20 ~ 1.64	mmol/L
	126	106 ~ 145	mg/dL
CRE	88	53 ~ 122	µmol/L
	0.99	0.60 ~ 1.38	mg/dL
TP	59	53 ~ 65	g/L
	5.9	5.3 ~ 6.5	g/dL
ALB	36	29 ~ 42	g/L
	3.6	2.9 ~ 4.2	g/dL

TBIL	14	10 ~ 19	μmol/L
	0.8	0.6 ~ 1.1	mg/dL
HDL-C	1.09	0.88 ~ 1.27	mmol/L
	42	34 ~ 49	mg/dL
Ca	2.15	1.83 ~ 2.48	mmol/L
	8.6	7.3 ~ 9.9	mg/dL
IP	1.16	1.00 ~ 1.29	mmol/L
	3.6	3.1 ~ 4.0	mg/dL
DBIL	6.8	3.4 ~ 12.0	μmol/L
	0.4	0.2 ~ 0.7	mg/dL
Mg	0.66	0.53 ~ 0.78	mmol/L
	1.6	1.3 ~ 1.9	mg/dL
GGT	63	53 ~ 74	U/L
GOT/AST	35	28 ~ 41	U/L
GPT/ALT	39	33 ~ 46	U/L
CPK	70	54 ~ 86	U/L
LDH	167	139 ~ 194	U/L
ALP	198	166 ~ 230	U/L
AMYL	61	49 ~ 73	U/L
LAP	50	35 ~ 66	U/L
CHE	283	245 ~ 320	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



**FUJIFILM Corporation**

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN