

QP-L :Mean & expected ranges

LOT 21355



2026.08

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	6.0	5.3 ~ 6.7	mmol/L
	108	95 ~ 120	mg/dL
BUN	5.00	4.32 ~ 5.68	mmol/L
	14.0	12.1 ~ 15.9	mg/dL
UA	220	184 ~ 250	µmol/L
	3.7	3.1 ~ 4.2	mg/dL
TCHO	3.28	2.90 ~ 3.65	mmol/L
	127	112 ~ 141	mg/dL
TG	1.49	1.26 ~ 1.72	mmol/L
	132	112 ~ 152	mg/dL
CRE	86	51 ~ 120	µmol/L
	0.97	0.58 ~ 1.36	mg/dL
TP	58	52 ~ 64	g/L
	5.8	5.2 ~ 6.4	g/dL
ALB	36	29 ~ 42	g/L
	3.6	2.9 ~ 4.2	g/dL

TBIL	12	7 ~ 17	μmol/L
	0.7	0.4 ~ 1.0	mg/dL
HDL-C	0.93	0.67 ~ 1.16	mmol/L
	36	26 ~ 45	mg/dL
Ca	2.35	2.03 ~ 2.68	mmol/L
	9.4	8.1 ~ 10.7	mg/dL
IP	1.19	1.03 ~ 1.32	mmol/L
	3.7	3.2 ~ 4.1	mg/dL
DBIL	5.1	1.7 ~ 10.3	μmol/L
	0.3	0.1 ~ 0.6	mg/dL
Mg	0.62	0.49 ~ 0.70	mmol/L
	1.5	1.2 ~ 1.7	mg/dL
GGT	63	53 ~ 74	U/L
GOT/AST	37	30 ~ 43	U/L
GPT/ALT	41	35 ~ 47	U/L
CPK	73	57 ~ 89	U/L
LDH	154	128 ~ 180	U/L
ALP	207	172 ~ 242	U/L
AMYL	54	42 ~ 67	U/L
LAP	50	35 ~ 65	U/L
CHE	271	235 ~ 307	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN