

QP-H :Mean & expected ranges

LOT 22356


2027.02

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.9	16.2 ~ 19.6	mmol/L
	322	291 ~ 353	mg/dL
BUN	17.92	16.14 ~ 19.67	mmol/L
	50.2	45.2 ~ 55.1	mg/dL
UA	648	595 ~ 702	µmol/L
	10.9	10.0 ~ 11.8	mg/dL
TCHO	3.85	3.44 ~ 4.29	mmol/L
	149	133 ~ 166	mg/dL
TG	2.63	2.35 ~ 2.91	mmol/L
	233	208 ~ 258	mg/dL
CRE	454	402 ~ 507	µmol/L
	5.14	4.55 ~ 5.73	mg/dL
TP	74	68 ~ 80	g/L
	7.4	6.8 ~ 8.0	g/dL
ALB	47	41 ~ 53	g/L
	4.7	4.1 ~ 5.3	g/dL

TBIL	75	62 ~ 89	μmol/L
	4.4	3.6 ~ 5.2	mg/dL
HDL-C

Ca	3.38	2.98 ~ 3.75	mmol/L
	13.5	11.9 ~ 15.0	mg/dL
IP	2.07	1.84 ~ 2.26	mmol/L
	6.4	5.7 ~ 7.0	mg/dL
DBIL

Mg

GGT	704	600 ~ 808	U/L
GOT/AST	101	90 ~ 112	U/L
GPT/ALT	109	96 ~ 121	U/L
CPK	279	241 ~ 316	U/L
LDH	395	341 ~ 450	U/L
ALP	650	562 ~ 738	U/L
AMYL	178	152 ~ 203	U/L
LAP	82	67 ~ 97	U/L
CHE	343	297 ~ 388	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları" nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN