

QP-H :Mean & expected ranges

LOT 22457


2027.08

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	17.6	15.9 ~ 19.3	mmol/L
	317	287 ~ 348	mg/dL
BUN	17.89	16.14 ~ 19.67	mmol/L
	50.1	45.2 ~ 55.1	mg/dL
UA	648	595 ~ 702	µmol/L
	10.9	10.0 ~ 11.8	mg/dL
TCHO	4.14	3.70 ~ 4.60	mmol/L
	160	143 ~ 178	mg/dL
TG	2.59	2.30 ~ 2.87	mmol/L
	229	204 ~ 254	mg/dL
CRE	451	399 ~ 503	µmol/L
	5.10	4.51 ~ 5.69	mg/dL
TP	74	67 ~ 80	g/L
	7.4	6.7 ~ 8.0	g/dL
ALB	42	35 ~ 48	g/L
	4.2	3.5 ~ 4.8	g/dL

TBIL	96	82 ~ 109	μmol/L
	5.6	4.8 ~ 6.4	mg/dL
HDL-C

Ca	3.20	2.83 ~ 3.60	mmol/L
	12.8	11.3 ~ 14.4	mg/dL
IP	2.03	1.84 ~ 2.26	mmol/L
	6.3	5.7 ~ 7.0	mg/dL
DBIL

Mg

GGT	680	580 ~ 781	U/L
GOT/AST	110	98 ~ 122	U/L
GPT/ALT	120	106 ~ 134	U/L
CPK	311	269 ~ 353	U/L
LDH	404	349 ~ 459	U/L
ALP	639	552 ~ 726	U/L
AMYL	172	147 ~ 197	U/L
LAP	88	73 ~ 103	U/L
CHE	360	312 ~ 407	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN