

## QP-L :Mean & expected ranges

**LOT** 21356


2027.02

Applicable to FDC analyzer using slides.  
 Applicable à l'analyseur FDC utilisant des lames.  
 Applicabile ad analizzatori FDC mediante slide.  
 Anwendbar auf den FDC-Analysator mit Objektträgern.  
 Aplicable al analizador FDC con slides.  
 Aplicável ao analisador FDC utilizando lâminas.  
 Применимо к анализатору FDC на основе слайдов.  
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	5.9	5.3 ~ 6.7	mmol/L
	107	95 ~ 120	mg/dL
BUN	4.96	4.32 ~ 5.64	mmol/L
	13.9	12.1 ~ 15.8	mg/dL
UA	220	184 ~ 250	µmol/L
	3.7	3.1 ~ 4.2	mg/dL
TCHO	3.13	2.77 ~ 3.49	mmol/L
	121	107 ~ 135	mg/dL
TG	1.50	1.28 ~ 1.73	mmol/L
	133	113 ~ 153	mg/dL
CRE	86	51 ~ 120	µmol/L
	0.97	0.58 ~ 1.36	mg/dL
TP	59	52 ~ 65	g/L
	5.9	5.2 ~ 6.5	g/dL
ALB	36	29 ~ 42	g/L
	3.6	2.9 ~ 4.2	g/dL

TBIL	14	9 ~ 19	μmol/L
	0.8	0.5 ~ 1.1	mg/dL
HDL-C	0.88	0.65 ~ 1.14	mmol/L
	34	25 ~ 44	mg/dL
Ca	2.35	2.03 ~ 2.68	mmol/L
	9.4	8.1 ~ 10.7	mg/dL
IP	1.16	1.03 ~ 1.29	mmol/L
	3.6	3.2 ~ 4.0	mg/dL
DBIL	6.8	3.4 ~ 12.0	μmol/L
	0.4	0.2 ~ 0.7	mg/dL
Mg	0.58	0.45 ~ 0.70	mmol/L
	1.4	1.1 ~ 1.7	mg/dL
GGT	61	51 ~ 72	U/L
GOT/AST	33	27 ~ 39	U/L
GPT/ALT	40	34 ~ 47	U/L
CPK	68	51 ~ 84	U/L
LDH	156	130 ~ 183	U/L
ALP	197	165 ~ 229	U/L
AMYL	54	42 ~ 66	U/L
LAP	50	35 ~ 65	U/L
CHE	273	237 ~ 309	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



**FUJIFILM Corporation**

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN