

## QP-L :Mean & expected ranges

**LOT** 21457


2027.08

Applicable to FDC analyzer using slides.  
 Applicable à l'analyseur FDC utilisant des lames.  
 Applicabile ad analizzatori FDC mediante slide.  
 Anwendbar auf den FDC-Analysator mit Objektträgern.  
 Aplicable al analizador FDC con slides.  
 Aplicável ao analisador FDC utilizando lâminas.  
 Применимо к анализатору FDC на основе слайдов.  
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campioni Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	5.7	5.1 ~ 6.4	mmol/L
	103	91 ~ 115	mg/dL
BUN	5.03	4.39 ~ 5.71	mmol/L
	14.1	12.3 ~ 16.0	mg/dL
UA	220	190 ~ 250	µmol/L
	3.7	3.2 ~ 4.2	mg/dL
TCHO	3.31	2.95 ~ 3.70	mmol/L
	128	114 ~ 143	mg/dL
TG	1.46	1.23 ~ 1.67	mmol/L
	129	109 ~ 148	mg/dL
CRE	88	53 ~ 122	µmol/L
	0.99	0.60 ~ 1.38	mg/dL
TP	58	52 ~ 64	g/L
	5.8	5.2 ~ 6.4	g/dL
ALB	31	25 ~ 37	g/L
	3.1	2.5 ~ 3.7	g/dL

TBIL	12	7 ~ 17	μmol/L
	0.7	0.4 ~ 1.0	mg/dL
HDL-C	0.93	0.67 ~ 1.16	mmol/L
	36	26 ~ 45	mg/dL
Ca	2.25	1.90 ~ 2.58	mmol/L
	9.0	7.6 ~ 10.3	mg/dL
IP	1.16	1.00 ~ 1.29	mmol/L
	3.6	3.1 ~ 4.0	mg/dL
DBIL	6.8	1.7 ~ 10.3	μmol/L
	0.4	0.1 ~ 0.6	mg/dL
Mg	0.58	0.45 ~ 0.70	mmol/L
	1.4	1.1 ~ 1.7	mg/dL
GGT	61	50 ~ 71	U/L
GOT/AST	39	33 ~ 46	U/L
GPT/ALT	46	39 ~ 52	U/L
CPK	73	56 ~ 89	U/L
LDH	162	135 ~ 189	U/L
ALP	185	153 ~ 217	U/L
AMYL	53	41 ~ 66	U/L
LAP	53	38 ~ 68	U/L
CHE	286	248 ~ 324	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



<https://med-lib.fujifilm.com/>



**FUJIFILM Corporation**

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN