

QP-L :Mean & expected ranges

LOT 21560


2028.07

Applicable to FDC analyzer using slides.
 Applicable à l'analyseur FDC utilisant des lames.
 Applicabile ad analizzatori FDC mediante slide.
 Anwendbar auf den FDC-Analysator mit Objektträgern.
 Aplicable al analizador FDC con slides.
 Aplicável ao analisador FDC utilizando lâminas.
 Применимо к анализатору FDC на основе слайдов.
 Kızaklar kullanan FDC analizör için geçerlidir.

Test item Élément de test Campione Testelement Elemento de prueba Item de teste Объект испытания Test ögesi	Mean Moyenne Valore medio Mittelwert Media Média Среднее Ortalama	Expected range Plage attendue Intervallo previsto Erwarteter Bereich Rango esperado Intervalo esperado Ожидаемые пределы Beklenen aralık	Unit Unité Unità Einheit Unidad Unidade Единица Birim
GLU	5.9	5.2 ~ 6.6	mmol/L
	106	94 ~ 118	mg/dL
BUN	5.03	4.39 ~ 5.71	mmol/L
	14.1	12.3 ~ 16.0	mg/dL
UA	214	184 ~ 244	µmol/L
	3.6	3.1 ~ 4.1	mg/dL
TCHO	3.26	2.90 ~ 3.62	mmol/L
	126	112 ~ 140	mg/dL
TG	1.47	1.24 ~ 1.68	mmol/L
	130	110 ~ 149	mg/dL
CRE	84	50 ~ 118	µmol/L
	0.95	0.56 ~ 1.34	mg/dL
TP	58	52 ~ 64	g/L
	5.8	5.2 ~ 6.4	g/dL
ALB	33	27 ~ 39	g/L
	3.3	2.7 ~ 3.9	g/dL

TBIL	12	7 ~ 17	μmol/L
	0.7	0.4 ~ 1.0	mg/dL
HDL-C	0.98	0.72 ~ 1.22	mmol/L
	38	28 ~ 47	mg/dL
Ca	2.20	1.88 ~ 2.53	mmol/L
	8.8	7.5 ~ 10.1	mg/dL
IP	1.16	1.03 ~ 1.32	mmol/L
	3.6	3.2 ~ 4.1	mg/dL
DBIL	6.8	1.7 ~ 10.3	μmol/L
	0.4	0.1 ~ 0.6	mg/dL
Mg	0.58	0.45 ~ 0.70	mmol/L
	1.4	1.1 ~ 1.7	mg/dL
GGT	63	52 ~ 73	U/L
GOT/AST	35	29 ~ 41	U/L
GPT/ALT	43	37 ~ 50	U/L
CPK	77	61 ~ 93	U/L
LDH	162	135 ~ 189	U/L
ALP	197	165 ~ 229	U/L
AMYL	58	46 ~ 70	U/L
CHE	276	240 ~ 313	U/L

Read "Instructions for use" before use.

Lire les « Instructions d'utilisation » avant utilisation.

Prima dell'uso leggere "Istruzioni per l'uso".

Vor der Verwendung Gebrauchsanweisung beachten.

Leer "Instrucciones de uso" antes del uso.

Leia as "Instruções de utilização" antes de utilizar.

Перед использованием прочитайте "Инструкции по применению".

Kullanmadan önce "Kullanım talimatları"nı okuyun.



FUJIFILM Corporation

26-30, Nishiazabu, 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN